

## Mouse Cell Line

		1000	1001
	Agents	Premium	Regular
ECTR	Ectromelia Virus	•	•
HAN	Hantaan Virus	•	
KV	K Virus (Mouse Pneumonitis Virus)	•	
LDHV	Lactate Dehydrogenase Elevating Virus	•	•
LCMV	Lymphocytic Choriomeningitis Virus	•	•
MMV	Mouse Minute virus	•	•
MAV-1	Mouse Adenovirus FL	•	•
MAV-2	Mouse Adenovirus K87	•	•
MPV	Mouse Parvovirus	•	•
EDIM	Rotavirus/Epizootic Diarrhea of Infant Mice	•	•
MHV	Mouse Hepatitis Virus	•	•
MYCO	<i>Mycoplasma</i> generic	•	
MCMV	Murine Cytomegalovirus	•	
MTV	Mouse Thymic Virus	•	
PVM	Pneumonia Virus of Mice	•	•
POLY	Polyoma Virus	•	
REO3	Respiratory Enteric Virus III	•	•
SEN	Sendai Virus	•	•
TMEV/GD7	Theiler's Mouse Encephalomyelitis Virus	•	•

## Rat Cell Line

		1002
	Agents	
HAN	Hantaan Virus	•
KRV	Killham Rat Virus	•
PVM	Pneumonia Virus of Mice	•
RCV	Rat Coronavirus	•
REO3	Respiratory Enteric Virus III	•
SEN	Sendai Virus	•
TMEV/GD7	Theiler's Mouse Encephalomyelitis Virus	•
H-1	Toolan's H-1 Virus	•

## Hamster Cell Line

		1003
	Agents	
PVM	Pneumonia Virus of Mice	•
REO3	Respiratory Enteric Virus III	•
SEN	Sendai Virus	•
LCMV	Lymphocytic Choriomeningitis Virus	•

401 Professional Dr., Suite 210  
Gaithersburg, MD 20879  
P: 301-610-2521 | F: 240-686-6776

7540 Louis Pasteur, Suite 200  
San Antonio, TX 78229  
P: 877-615-7275 | F: 210-615-7771