

CytoID CELL LINE AUTHENTICATION (STR)

SAMPLE TYPE	TEST CODE	INCLUDES MYCOPLASMA SPP.
Human 9 STR marker profile: Interspecies contamination check for human, mouse, rat, African green monkey, and Chinese hamster	HSTR9	HSTM9
Human 16 STR marker profile: Interspecies contamination check for human, mouse, rat, African green monkey, and Chinese hamster	HST16	STM16
Mouse genetic profile using NIST published 9 species-specific STR markers: Interspecies contamination check for human, mouse, rat, African green monkey, and Chinese hamster	MSTR9	MSTM9
Rat genetic profile using NIST published 31 species-specific STR markers: Interspecies contamination check for human, mouse, rat, African green monkey, and Chinese hamster	STR31	STM31