

RODENT VIRUS/INDIVIDUAL PCR

ANIMAL HEALTH MONITORING

TEST CODE	ABBREVIATION	AGENT	SAMPLE TYPE	SPECIES
10569	BC	Boone Cardiovirus	Feces, environmental	Rat
10424	ECTRO	Ectromelia Virus	Feces, environmental, skin lesion	Mouse
10526	GPADV	Guinea pig Adenovirus	Feces, environmental, oral swab, lung	Guinea Pig
10425	GPCMV	Guinea pig cytomegalovirus	Feces, spleen, salivary gland	Guinea Pig
10426	HANT	Hantaan Virus	Feces, environmental, kidney	Mouse, rat, guinea pig, hamster
10401	PUUV	Puumala Virus	Feces, environmental, kidney	Mouse, rat
10402	SINV	Sin Nombre Virus	Feces, environmental, kidney	Mouse
10403	SEOV	Seoul Virus	Feces, environmental, kidney	Mouse, rat
10445	H-1	Toolan's H-1 Virus	Feces, environmental, mesenteric lymph nodes, spleen	Rat
10446	KRV	Killham Rat Virus	Feces, environmental, mesenteric lymph nodes, spleen	Rat
10449	KRV/H-1	Killham Rat Virus/Toolan's H-1 Virus	Feces, environmental, mesenteric lymph nodes, spleen	Rat
10433	KV	K Virus (Mouse Pneumonitis Virus)	Feces, environmental, mesenteric lymph nodes	Mouse
10434	LDHV	Lactate Dehydrogenase Elevating Virus	Serum, spleen	Mouse
10435	LCMV	Lymphocytic Choriomeningitis Virus	Feces, environmental, kidney	Mouse
10436	MAV-1	Mouse Adenovirus 1 (FL)	Feces, environmental, mesenteric lymph nodes	Mouse
10529	MAV-2	Mouse Adenovirus 2 (K87)	Feces, environmental, intestine	Mouse
10438	MHV	Mouse Hepatitis Virus	Feces, environmental, mesenteric lymph nodes	Mouse
10437	MCMV	Murine Cytomegalovirus	Feces, environmental, spleen, salivary gland	Mouse
10447	MVM (MMV)	Mouse Minute virus	Feces, environmental, mesenteric lymph nodes, spleen	Mouse
10448	MPV	Mouse Parvovirus	Feces, environmental, mesenteric lymph nodes, spleen	Mouse
10450	MMV/MPV	Mouse Minute virus/Mouse Parvovirus	Feces, environmental, mesenteric lymph nodes, spleen	Mouse
10439	EDIM	Rotavirus/Epizootic Diarrhea of Infant Mice	Feces, environment	Mouse
10460	MTV	Mouse Thymic Virus	Feces, environmental, oral swab, salivary gland	Mouse
10440	MNV	Murine Norovirus	Feces, environmental	Mouse
10444	PARVO	Parvovirus generic	Feces, environmental, mesenteric lymph nodes, spleen	Mouse, rat
10452	PVM	Pneumonia Virus of Mice	Feces, environmental, oral swab, lung	Mouse
10453	POLY	Polyoma Virus	Feces, kidney, environmental	Mouse
10454	POLY/K	Polyoma Virus / K virus	Feces, kidney, environmental	Mouse
10456	REO	Respiratory Enteric Virus III	Feces, liver, lung, environmental	Mouse, rat, hamster, guinea pig
10455	RCV	Rat Coronavirus	Feces, environmental, oral swab, Harderian gland	Rat
10587	RCMV	Rat Cytomegalovirus	Feces, environmental, Oral swab, salivary gland	Rat
10473	RMV	Rat Minute Virus	Feces, environmental, mesenteric lymph nodes, spleen	Rat
10474	RPV	Rat Parvovirus	Feces, environmental, mesenteric lymph nodes, spleen	Rat
10457	SEN	Sendai Virus	Feces, environmental, oral swab, lung	Mouse, rat, guinea pig, hamster
10459	TMEV/GDVII	Theiler's Mouse Encephalomyelitis Virus	Feces, environmental	Mouse

RODENT PARASITES/INDIVIDUAL PCR

TEST CODE	ABBREVIATION	AGENT	SAMPLE TYPE	SPECIES
10479	FUR	Fur mites	Fur swab, skin swab	Mouse, rat
10461	PIN	Pinworm (Syphacia & Aspicularis)	Feces, environmental	Mouse, rat
10478	CRYP	Cryptosporidium spp.	Feces, environmental	Mouse, rat
10494	EIME	Eimeria spp.	Feces, environmental	Rabbit
10493	EMUR	Entamoeba muris	Fecal, environmental	Mouse, rat
10520	RNAN	Rodentolepis nana	Fecal, environmental	Mouse, rat
10521	HDIM	Hymenolepis diminuta (rat tapeworm)	Fecal, environmental	Mouse, rat
10492	TMUR	Tritrichomonas muris	Fecal, environmental	Mouse, rat
10498	GMUR	Giardia muris	Fecal, environmental	Mouse, rat
10487	GIAR	Giardia spp.	Feces, environmental	Mouse, rat
10486	SMUR	Spironucleus muris	Fecal, environmental	Mouse, rat
10524	DEMO	Demodex	Fur swab, environmental, feces	Mouse, rat

RODENT BACTERIA/INDIVIDUAL PCR

TEST CODE	ABBREVIATION	AGENT	SAMPLE TYPE	SPECIES
10480	BBRO	Bordetella bronchiseptica	Feces, environmental, oral swab, trachea	Mouse, guinea pig, rat, rabbit
10404	BHIN	Bordetella hinzii	Feces, environmental, oral swab, trachea	Mouse
10405	CAMPY	Campylobacter spp	Feces, environmental	Mouse, rat, hamster, guinea pig
10406	CCOL	Campylobacter coli	Feces, environmental	Mouse, rat, hamster, guinea pig
10407	CJEJ	Campylobacter jejuni	Feces, environmental	Mouse, rat, hamster, guinea pig
10424	CARB	Cilia Associated Respiratory Bacillus	Feces, environmental, oral swab, trachea	Mouse, rat, rabbit
10408	CROD	Citrobacter rodentium	Feces, environmental	Mouse, rat
10590	CPER	Clostridium perfringens	Feces, environmental	Mouse, rat
10421	CPILI	Clostridium piliforme	Fecal, environment	Mouse, rat, rabbit, ferret
10422	CBOV	Corynebacterium bovis	Feces, environmental, fur swab	Mouse, rat
10409	CKUT	Corynebacterium kutscheri	Feces, environment, oral swab	Mouse, hamster, rat
10420	ECUN	Encephalitozoon cuniculi	Feces, kidney, brain, environmental	Mouse, rat, hamster, guinea pig, rabbit
10427	HELICO	Helicobacter generic	Feces, environmental	Mouse, guinea pig, rabbit, rat
10428	HBIL	Helicobacter bilis	Feces, environmental	Mouse, rabbit, rat
10429	HHEP	Helicobacter hepaticus	Feces, environmental	Mouse, rabbit, rat
10462	HGAN	Helicobacter ganmani	Feces, environment	Mouse, rabbit, rat
10525	HMASTO	Helicobacter mastomyrinus	Feces, environment	Mouse, rabbit, rat
10430	HROD	Helicobacter rodentium	Feces, environment	Mouse, rabbit, rat
10431	HTYP	Helicobacter typhlonicus	Feces, environment	Mouse, rabbit, rat
10432	HMUL	Helicobacter Multiplex	Feces, environmental	Mouse, rabbit, rat
10410	LINT	Lawsonia intracellularis	Feces, intestine, mesenteric lymph node	guinea pig, hamster, rat, rabbit
10573	LEPTO	Leptospira	Feces, kidney	Mouse, rat
10411	KOXY	Klebsiella oxytoca	Feces, environment	Mouse, rat
10412	KPNE	Klebsiella pneumoniae	Feces, environment	Mouse, rat
10441	MycO spp	Mycoplasma spp	Feces, environmental, oral, swab	Mouse, rat
10442	MART	Mycoplasma arthritidis	Feces, environmental, oral, swab	Mouse, rat
10443	MPUL	Mycoplasma pulmonis	Feces, environmental, oral, swab	Mouse, rat
14013	PMUL	Pasteurella multocida	Feces, environmental, oral, swab	Rat, mouse, rabbit, guinea pig, hamster
10531	PAST	Pasteurella spp.	Feces, environmental, oral, swab	Rat, mouse, rabbit, guinea pig, hamster
10451	PCAR	Pneumocystis carinii	Oral swab, lung, environmental	Rat
10414	PMUR	Pneumocystis murina	Oral swab, lung, environmental	Mouse
10415	PNEU	Pneumocystis spp	Oral swab, lung, environmental	Rat, mouse
10416	PMIR	Proteus mirabilis	Feces, environment	Rat, mouse, rabbit, guinea pig, hamster
10417	PAER	Pseudomonas aeruginosa	Feces, environmental	Mouse, rat
10463	RPNEU	Rodentibacter pneumotropica	Feces, environmental, oral swab, trachea	Mouse, rat
10579	RHEYL	Rodentibacter heylII	Feces, environmental, oral swab, trachea	Mouse, rat
10418	SALM	Salmonella spp.	Feces, environmental	Mouser, guinea pig, hamster, rabbit, rat
10588	SFB	Segmented Filamentous Bacteria	Feces, environmental	Mouse
10458	SAUR	Staphylococcus aureus	Feces, environmental, oral swabs, fur swab/skin swab	Mouse, rat
10522	STAX	Staphylococcus xylosus	Feces, environmental, oral swab, fur swab/skin swab	Mouse, rat
10464	SMONI	Streptobacillus moniliformis	Feces, environmental, oral swab	Mouse, rat
10523	SBETA	B-hemolytic Streptococcus (Groups A, B, C, G)	Feces, environmental, oral swab	Mouse, rat
10419	SAGA	Streptococcus agalactiae	Feces, environmental, oral swab	Mouse, rat
10495	SCAN	Streptococcus canis	Feces, environmental, oral swab	Mouse, rat
10496	SEQU	Streptococcus equi subs. Zooepidemicus	Feces, environmental, oral swab	Mouse, rat
10465	SPNE	Streptococcus pneumoniae	Feces, environmental, oral swab	Mouse, rat
10497	SPYO	Streptococcus pyogenes	Feces, environmental, oral swab	Mouse, rat
10482	TREP	Treponema cuniculi	Genital swab	Rabbit